


iCS-digital™ PSC kit Troubleshooting Guide

I. QuantaSoft™ Analysis Pro results troubleshooting table

Problem observed		Possible reason	Solution	
<p>Small 2D amplitude clusters and low concentrations (<150 copies/μL)</p> <p><i>See Fig.1. as an example</i></p>		<ul style="list-style-type: none"> - Low DNA concentrations were used 	<ul style="list-style-type: none"> - Check the concentration of your DNA using a Qubit fluorometer and ensure you dilute the DNA samples to a concentration of 5 ng/μL 	
<p>Large 2D amplitude clusters & high concentrations (>250 copies/μL)</p> <p><i>See Fig.2. as an example</i></p>		<ul style="list-style-type: none"> - High DNA concentrations were used 		
<p>Intense Rain and/or duplicated clusters make the analysis difficult or impossible</p> <p><i>See Fig.3. as an example</i></p>	<p>All wells affected</p>		<ul style="list-style-type: none"> - The Mix assays were not vortexed enough - PCR program is incorrect 	<ul style="list-style-type: none"> - Vortex vigorously each Mix assay tube twice for 5-10 seconds, and briefly centrifuge the tubes between each vortexing steps - Check that the PCR program conforms to the protocol
	<p>Only few wells affected</p>	<p>droplet number < 10 000</p>	<ul style="list-style-type: none"> - Probable issue with the QX instrument (dust, bubbles, etc.) 	<ul style="list-style-type: none"> - Perform the test in duplicate
		<p>droplet number ≥ 10 000</p>	<ul style="list-style-type: none"> - Some Mix assays were not vortexed enough 	<ul style="list-style-type: none"> - Vortex vigorously each Mix assay tube twice for 5-10 seconds, and briefly centrifuge the tubes between each vortexing steps


Fig.1. Small 2D amplitude clusters and/or low concentrations (<200 copies/ μ L)



A. 2D Amplitude tab. B. Concentration tab (normal values \sim 200 copies/ μ L).

Left panels: 2D amplitude & Concentration tab of a sample with correct DNA concentration. Right panels: 2D amplitude & Concentration tab of a sample with low DNA concentration.


Fig.2. Large 2D amplitude clusters and/or high concentrations (>200 copies/ μ L)



A. 2D Amplitude tab. B. Concentration tab (normal values ~200 copies/ μ L).

Left panels: 2D amplitude & Concentration tab of a sample with correct DNA concentration. Right panels: 2D amplitude & Concentration tab of a sample with high DNA concentration.

Fig.3. Intense rain and/or duplicated clusters make the analysis difficult or impossible



Example of a low quality 2D-cluster well

II. Online iCS-digital™ analysis tool troubleshooting table

Problem observed	Possible reason	Problem resolution
<p>“502 bad gateway nginx/1.19.6” Error message</p>	<ul style="list-style-type: none"> - Incorrect naming of the Target Example: “3q” instead of “3p” 	<ul style="list-style-type: none"> - Validate the correct naming of the targets (column C of the Input Excel file) Refer to Page 6 & Pages 11-12 (APPENDIX 2) of the User Manual
	<ul style="list-style-type: none"> - Use of a special character for the Sample name (e.g., ~ ! @ # \$ ^ % & * ? { }) 	<ul style="list-style-type: none"> - Remove any special characters present in the Sample names (column B of the Input Excel file) Refer to Page 6 of the User Manual
	<ul style="list-style-type: none"> - Incorrect sorting of the data prior to their exportation from the QuantaSoft™ Analysis Pro software 	<ul style="list-style-type: none"> - In the “2D Amplitude” tab of the QuantaSoft™ Analysis Pro software, sort the lines by “Sample” name in the Well Data table prior to exporting the data in Excel format Refer to Page 9 (Figure 3) & Pages 13-14 (APPENDIX 3) of the User Manual
	<ul style="list-style-type: none"> - A sample with concentration (copies/μL) and therefore CNV values equal to zero (e.g., used of a NCT) 	<ul style="list-style-type: none"> - Do not include samples with concentration (column “Conc” of the Well Data table) equal or close to zero when exporting the data from the QuantaSoft™ Analysis Pro software
<p>Incoherent sample name in the final report</p>	<ul style="list-style-type: none"> - Incorrect sorting of the data prior to their exportation from the QuantaSoft™ Analysis Pro software 	<ul style="list-style-type: none"> - In the “2D Amplitude” tab of the QuantaSoft™ Analysis Pro software, sort the lines by “Sample” name in the Well Data table prior to exporting the data in Excel format Refer to Page 9 (Figure 3) & Page 13-14 (APPENDIX 3) of the User Manual